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1.1 Thank you very much

Thank you for your confidence in us and congratulations!

The **PSC** personal scales with body analysis functionality which you have purchased are a top quality product from **MEDISANA**.

This unit is intended for weighing and calculating the proportion of body fat, water content and muscle mass in people.

In order to achieve the desired effect with your **MEDISANA** Personal scales **PSC** with body analysis functionality in the long term, we recommend that you read the following information on its use and maintenance carefully.

1.2 Important

Please read the following safety information carefully before using the device and keep these instructions in a safe place for later use.

- Any treatment or diet for people who are underweight or overweight requires the qualified advice of a professional (doctor, dietician). The values determined with the scales may be source of useful support.
- This device is not suitable for pregnant women!
- These scales are not suitable for people with pacemakers or other medical implants. The body fat information may be inaccurate in people with diabetes or other medical/physical limitations. The same also applies to very highly trained athletes.

- ***Risk of tipping over!***

Do not place the scale on an uneven surface.

Do not stand on the corners/edges of the scale.

When weighing yourself, do not stand on one side or the edge of the scales.

- ***Risk of slipping!***

Do not stand on the scale with wet feet.

Do not stand on the scale wearing socks.

- Use the device only according to its intended purpose as specified in the instruction manual.
- The warranty will be invalidated if the device is used for purposes other than those for which it is intended.
- The scale has been manufactured for domestic use. It is not suitable for commercial use in hospitals or other medical institutions.
- Do not use the device if it is not working correctly, if it has been dropped or has fallen into water or has been damaged.
- The scales have a measuring range up to 180 kg, 400 lbs. Do not overload the scales.
- Place the scales on a firm, level base. Soft, uneven surfaces are unsuitable for the measuring process and lead to incorrect results.
- Do not use the scale in a location where extremely high temperatures or humidity may occur.
- Keep the scale away from water.

1.2 Important

- Handle the scale with care. The scale must not be subjected to impacts or vibration. Do not drop the scale.
- Do not attempt to dismantle the scale, otherwise the warranty will be invalidated. The scale does not contain any parts that need to be maintained by the user, with the exception of the battery.
- Do not attempt to repair the device yourself in the event of a malfunction since this will invalidate the warranty. Repairs should only be carried out by authorised service centres.
- If the scale is not going to be used for a long period, please remove the batteries.

1.3 Scope of supply and packaging

Please check first of all that the unit is complete.

The following parts are included as standard:

- 1 **MEDISANA PSC** Personal scales with body analysis functionality
- 2 Type CR2032 3V lithium coin cell batteries
- 1 Instructions for use

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

WARNING

To avoid risk of suffocation, keep packaging film away from children!

2.1 Performance features

- Extremely flat design, only 20 mm high
- Top quality safety glass
- Top quality sensor technology
- Age entry from 10 – 99 years
- Height entry from 75 to 215 cm
- Measuring range 2 kg to 180 kg or 4.4 lbs to 400 lbs
- 0.1 kg or 0.2 lb graduation
- Switchable between KG, LB and ST
- Body fat analysis in 0.1 % steps
- Water content measurement in 0.1 % steps
- Muscle mass measurement in 0.1 % steps
- 12 memory locations
- Direct measurement of body weight
- Easy to use
- Clear LCD display
- Includes 2 type CR2032 3V batteries

2.2 Why is it important to know how much body fat you have?

The human body is made up amongst other things of water, fat and muscle mass. A fat proportion which is too high or too low may lead to increased health risks. The personal scales with body analysis functionality provide useful support in motivating someone to reach or maintain the best possible body weight by monitoring the proportion of body fat.

2.3 How does the scale with body analysis functionality work?

The scale measures using BIA (Bio-electrical Impedance Analysis). When you step onto the metal plates barefoot, extremely weak electrical current is passed into your body. The scale weighs the difference between these signals. The measurements are used in combination with personal data such as age, gender, height and weight to calculate the amount of body fat.

2.4 Prerequisites for correct measuring results

NOTES

Excessive drinking, eating, training, medical treatment, the female menstrual cycle etc. affect the measuring results. The body fat, water content and muscle measurements will be incorrect for women during pregnancy.

Reliable measurements can only be achieved under certain prerequisites:

- Always measure barefoot.
- The feet must have a good contact with the electrodes on the standing area.
- Stand upright.
- Carry out the measurement approximately 15 minutes after getting up, or after bathing or showering, with towel-dried feet.
- Please measure at the same time of day if possible.
- Do not measure after strenuous physical exercise. Wait at least an hour before measuring.
- There must not be any dirt stuck to the soles of your feet. Erroneous measuring results may occur if the soles of your feet are dry or extremely callused.

The amount of body fat, water content and muscle mass can only be reliably monitored using the measurements if they are taken consistently. Measurements that are taken under varying conditions can return different values.

 **NOTES**

The information in the tables is for orientation purposes only. If you have any questions about your proportion of body fat, please consult your doctor.

**2.5
Body fat proportion
range in %**

<i>Female</i>	Age	under-weight	normal	over-weight	fat
	15-39	<20,0	20,0-32,4	32,5-38,0	>38,0
	40-55	<23,0	23,0-35,4	35,5-41,0	>41,0
	>56	<24,0	24,0-36,4	36,5-42,0	>42,0
<i>Male</i>	Age	under-weight	normal	over-weight	fat
	15-39	<12,0	12,0-20,4	20,5-25,0	>25,0
	40-55	<13,0	13,0-21,4	21,5-26,0	>26,0
	>56	<14,0	14,0-22,4	22,5-27,0	>27,0

Body fat value analysis:

In order to analyse your body fat proportion, please compare your value to the values in the body fat table.

**2.6
Body water
proportion range
in %**

<i>Female</i>	Age	low	normal	high
	10-15	< 57	57-67	> 67
	16-30	< 47	47-57	> 57
	31-60	< 42	42-52	> 52
	61-80	< 37	37-47	> 47
<i>Male</i>	Age	low	normal	high
	10-15	< 58	58-72	> 72
	16-30	< 53	53-67	> 67
	31-60	< 47	47-61	> 61
	61-80	< 42	42-56	> 56

Body water value analysis:

In order to analyse your body water value, please compare your value with the values in the body water table.

**2.7
Average
Muscle Mass**

The average proportion of muscle mass for women between 15 and 70 years of age is max. 34% and for men 40%. Compare your value with the standard value to analyse your muscle mass.

3.1 Insert/replace batteries

Before you start up your digital scale, insert the two provided cell batteries (type CR2032) into the device. To do this, open the battery compartment **2** beneath the device and insert the batteries. Please note polarity (positive pole upwards; facing towards the lid). Close the battery compartment again. If the batteries are already inserted, just pull out the insulating strip. Replace the batteries if the battery change symbol "Lo" appears in the display **1** or if nothing appears in the display when the device is switched on.

BATTERY SAFETY INFORMATION

- **Keep away from children.**
- **Do not short-circuit.**
- **Do not recharge.**
- **Do not throw into fire.**
- **Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.**

3.2 Weight only

1. Place the scale on a stable, firm surface. All four feet must be touching the floor.
2. Touch the Start button **7** to switch on the scales.
3. When "0.0 kg" appears in the display **1**, step on the scales and stand still.
4. Your weight is measured and displayed
5. Step down from the scales. The scales switch off automatically after 20 seconds if you do not do anything else.

3.3 Change the Unit of Weight

When the scales are switched off, you can use the switch underneath to set the units of weight (*KG - LB - ST*).

3.4 Programming personal data

For more accurate readings, you must enter your personal information before weighing yourself. The scales can store personal information for up to 12 people: sex, height and age. The information is retained until new data is entered.

1. Place the scale on a stable, firm surface. All four feet must be touching the floor.
2. Press the ON button **3** to enter the information. The preset number "01" appears in the display. Press the Up (+) button **4** or the DOWN (-) button **5** to select the preset number. Press the OK button **6** to confirm your selection. The display switches to the sex setting.

3. Enter your sex by pressing the Up (+) button ④ or the DOWN (-) button ⑤ to select the appropriate symbol in the display. Press the OK button ⑥ to confirm your selection. The display switches to the height setting.
4. Press the Up (+) button ④ or the DOWN (-) button ⑤ to enter your height. You can enter a height from 75 to 215 cm in 1 cm increments. Press the OK button ⑥ to confirm your selection. The display switches to the age setting.
5. Press the Up (+) button ④ or the DOWN (-) button ⑤ to enter your age. You can enter an age from 10 - 99 years. Press the OK button ⑥ to confirm your selection. The preset number "01" appears again in the display.
6. The setting procedure is now complete. The scales switch off automatically after 3 seconds.

3.5 Weighing and measuring body fat, water content and muscle mass

NOTE

Reliable measurements can only be achieved under certain circumstances, as described in "2.4 Prerequisites for correct measuring results", p. 12.

1. Place the scales on a firm, hard surface. All four feet must be in contact with the floor.
2. Touch the Start button ⑦ and keep it pressed down until your preset number appears in the display. The readings last stored for you appear: weight, body fat, body water and body muscle. Before the first measurement, "888" appears in the display. After 2 seconds, "0.0" appears in the display.
3. Step onto the scales with bare feet and stand still. Place your feet on the electrodes ⑧. First, your weight is displayed. A moving "□" then appears in the display, indicating that your body fat, water and muscle are now being automatically measured. When the "□" disappears from the display, your total body fat (TBF), body water (TBW) and body muscle (TBM) are displayed.

Displays

TBF (Total Body Fat) = measured body fat
fat
overweight
normal
slim

TBW (Total Body Water) = measured body water

TBM (Total Body Muscle) = measured body muscle

4. The measurement is complete. Step down from the scales.
5. The display **1** once again shows – in succession – the measured body fat, body water and body muscle. Then the scales switch off automatically.

3.6 Error messages

A rectangular digital display showing the text "O-Ld" in a black, pixelated font.

The scale is overloaded.

A rectangular digital display showing the text "Lo" in a black, pixelated font.

The battery is low and must be replaced.

A rectangular digital display showing the text "Err" in a black, pixelated font.

Error during measurement.

4.1 Cleaning and maintenance

- Remove the batteries before cleaning the unit.
- Never use aggressive cleaning agents or stiff brushes.
- Clean the scale using a soft cloth lightly moistened with a mild soapy solution. Never use abrasive cleaning agents or alcohol.
- Do not allow any water to get into the unit. Do not use the unit again until it is completely dry.
- Remove the batteries from the unit if it will be unused for a longish period. Otherwise there is a risk that the batteries might leak.
- Do not expose the unit to direct sunlight; protect it against dirt and moisture.

4.2 Disposal

This product must not be disposed together with the domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Please remove the batteries before disposing of the equipment. Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.

Consult your municipal authority or your dealer for information about disposal.

4.3 Specifications

Name and model	:	MEDISANA PSC Personal scales with body analysis functionality
Dimensions (L x W x H)	:	approx. 30 cm x 31 cm x 2 cm
Weight	:	approx. 2.1 kg
Voltage supply	:	6 V= , 2 x 3 V lithium coin cell batteries type CR2032
Display system	:	Digital display
Memory	:	for 12 persons
Measuring range	:	2 kg to 180 kg or 4.4 lbs to 400 lbs
Maximum measuring deviation	:	± 1 %
Body fat analysis (TBF)	:	in steps of 0.1 %
Body water (TBW)	:	in steps of 0.1 %
Muscle mass (TBM)	:	in steps of 0.1 %
Graduation	:	100 g or 0.2 lb
Autom. switch-off	:	after approx. 3 sec.
Item number	:	40430
EAN number	:	4015588 40430 6

In accordance with our policy of continual product improvement, we reserve the right to technical and design changes without notice.

5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
 - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
 - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

MEDISANA AG

Itterpark 7-9
D40724 Hilden

Tel.: +49 (0) 2103 / 2007-60
Fax: +49 (0) 2103 / 2007-626
eMail: info@medisana.de
Internet: www.medisana.de

The service centre address is shown on the last page.

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D&L Marketing Oy
ALERE Electronics AB